



Nano multi needle 9pins

34G 

Decrease and minimize the sensation relative to single needle and thicker needle



More sanitary treatment with disposable suction filter

Increase the speed of the procedure

Covers larger treatment area

Less leakage of drug

* 32G 9pins optional



new skin booster injector

skin booster injections restores volumes and enhances skins glow

BellaVita™ Aqua

Cross Linked Hyaluronic Acid 20mg with Lidocaine & Mannitol

100% CROSS-LINKED HYALURONIC ACID
<3syringes / Box>



AQUA VIAL
<5vials / Box>



BELLAVITA™
P D R N

• PDRN 2% FROM SALMON SPERM
• CROSS-LINKED HA 10mg/cc
<2syringes / Box>



PDRN·S
by BELLAVITA

• PDRN 0.5% FROM SALMON SPERM
• CROSS-LINKED HA 5mg/cc
<5vials / Box>



CONTACTAC
CONTACTACOREA Inc.
Creative Medical & Cosmetic Service

Manufactured by CONTACTACOREA Inc.
#B119, B120 BizCenter-dong, SKn Technopark, 124,
Sagimakgol-ro, Jungwon-gu, Seongnam-si, Gyeonggi-do,
Korea. 13207
www.contactaco.com contactacorea@contactaco.com

Head Office - Seoul
1st Fl., SJ Bldg., 46, Nonhyun-ro 151-gil,
Gangnam-gu, Seoul, Korea. 06036
Tel : +82-2-540-6440 Fax : +82-2-540-6441



영문 | ENGLISH

bellapro

Bellapro Injector covers a larger treatment area, consistently delivering more accurate dosing and injections calibrated to a more precise depth than ever before.



System Feature

Two Precise Swiss Motor Installed

Motor, a heart of Injector, is upgraded by two swiss made motors. Calibrated to a more precision depth adjustment and accurate dosing by them.

Minimizes Drug Loss

To minimize drug loss newly upgraded technical program is applied.

Virus Sterilization Functions Installed

UV light cleaner installed kills 99.99% of viruses including CoronaVirus and bacteria in few seconds.

Customizing Treatment

Programed customizable treatment options for face, neck & décolletage and scalp.


Unique Vacuum Control For Less Pain

By low noise and soft vacuum pressure design, it is more quite, less bleeding and less pain.

Manual And Automatic Control Programed



Automatic injection program and removed suction function for scalp treatment

Wide 10.2 Inch Touch Screen




- 1 Mode : Face / Neck / Scalp (Filler, Solution, Filler + Solution)
- 2 Control Panel
- 3 Sound ON / OFF
- 4 UV LED light ON / OFF
- 5 Reset
- 6 Injection Counts
- 7 Standby / Ready

Human Centered Design



Programed treatment indicated

Programed customizable treatment options for face, neck& décolletage,scalp and hand.



Virus Sterilization Functions

Sterilize the surface of the multi needle that comes into contact with the skin. *Ref.



*References

1. Manuela B(2020) Far-UVC light efficiently and safely inactivates airborne Human coronavirus. Columbia University Irving Medical Center, New York, New York, SCIENTIFIC REPORTS 10:10285
2. N.Kazuhiro. "Disinfection application of UVC-LED", JUVA Newsletter,2018 No.11,13
3. ICNIRP Note on use of UVC lamps to kill/inactivate the coronavirus (SARS-CoV-2) 15 May 2020